

DSC-200T , Table Top Centrifuge

[Categories](#) > [Centrifuge Machine](#) > [Small Centrifuge](#) > [Table Top Centrifuge](#)



Product Name:

Table Top Centrifuge

Description:

Table Top Centrifuges are widely used in many fields such as biology, medicine, pharmacology, laboratory, and other regions.

Table Top Centrifuges are small, compact, strong and easy and safe for operation.
The max. Speed of Table Top Centrifuge is 3000 rpm at 60 Hz.

Features:

The design is up-to-date, beautiful outlook and durable.
The speed and time of Table Top Centrifuge are adjustable.
ABS outer shell can resist collision.
Brushless, free maintenance.
Table Top Centrifuge could accommodate 7ml, 10ml, and 15ml centrifuge tubes without changing the rotor. (7ml tube adapters are optional accessory)
With safety switch, when cover is opened, the motor can be stopped.
Speed is linear and stable.
Auto-balancing system.

Specification:

Model	Max. speed (rpm)	Timer (min.)	Speed & R.C.F. display	Manual brake	Motor	Weight (kg) (N.W. / G.W.)	Size (mm) (WxDxH)
DSC-200T	3,000 at 60Hz 2,500 at 50Hz	0~30 min. mechanical timer	None	None	Brushless induction motor	3.5 / 4.5	280×280×255

Specification of rotor for DSC-200T

Model	Max. Speed	Rotor type	Capacity
AR-1504	3,000 rpm	45 degree angle	4 tubes x (10~15)ml
AR-1506	3,000rpm	45 degerr angle	6 tubes x (10~15)ml
AR-1508	3,000rpm	45 dergerr angle	8 tubes x (10~15)ml

Note: 1.5ml~2.0ml and 5ml~7ml tube adapters are optional accessories

NOTE ; MAM E2 POH UNG PRICE NAMIN POH SA MGA CLINICAL CERTRIFUGE.

DSC-156T	Table Top Centrifuge, 6 Placer	12,250.00
DSC-158T	Table Top Centrifuge, 8 Placer	13,200.00
DSC-1512T	Table Top Centrifuge, 12 Placer	28,000.00