

H-500

Urine Analyzer



DIRUI H-500 Urine Analyzer's outstanding performance proves itself to be a reliable product for medium laboratories

Product Features:

- © 514 samples /hour throughput
- © 2000 patient results data memo
- © Easy operation by automatic detecting ready-to-test strips
- © Test items: Urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocytes, glucose, specific gravity, pH,VC, microalbumin,creatinine,calcium
- © Standard barcode reader
- ©Large LCD touch screen
- © With waste urine collecting system
- ©Built in thermal printer

Order Information:

Product Number	Description	Pack Size
1520104001	DIRUI H-500 Urine Analyzer	1 instrument

Product Number	Description	Volume	Pack Size
2310201003	DIRUI H8	100 strips/ bottle	1 bottle
2310201002	DIRUI H10	100 strips/ bottle	1 bottle
2310201001	DIRUI H11	100 strips/ bottle	1 bottle
2310201021	DIRUI H11MA	100 strips/ bottle	1 bottle
2310201037	DIRUI H13-Cr	100 strips/ bottle	1 bottle
2310201089	DIRUI H14-Ca	100 strips/ bottle	1 bottle